

PLANET LPI

ENCL



# UV LASER PRINTING

Exquisite and Fine Printing on Tablet and Capsule

**ENCLONY**


# FINE AND INNOVATION DESIGN



Fine printing technology allows you to print your ideal design. Lines can be printed with a minimum of 100µm, and the maximum surface of the tablet can be used for the design, such as the edges and corners.


## UV LASER PRINTING WITH FULL VISION INSPECTION

PROCESS




The machine performs UV laser printing with shape tracing first, and then operates full inspection. For inspection, each tablet is inspected by ten 2D cameras capturing the complete surface areas without any blind spot. In addition two 3D cameras measure the height of tablet surface, both sides of top and bottom, which allow to detect minimum 30µm defects.

## SHAPE TRACING AND KEY TECHNOLOGIES



- 1) Shape tracing
- 2) Laser Printing as per its shape
- 3) Inspection




In addition to precise shape tracing, a complex UV laser control method, and precise motor movement control, accurate printing locations are made within 100µm of the XY axis and 1° of rotation. Due to its ability to read the convex surface of the tablet well, the fine printing can be performed even on deep convex tablets without distortion.



# APPLICABLE FLEXIBILITY OF TABLET / CAPSULE

- Round tablet : Diameter(D) 3.5mm – 14mm
- Oval/Oblong tablet : Width(W) 3.5mm – 14mm  
Length(L) 10mm – 22mm
- Thickness(T) : 2mm – 8mm \*Diameter (or Width) minus Thickness  $\geq 0.7$ mm
- Polygon tablet : triangle, square, diamond, pentagon, hexagon, octagon shapes, etc.
- Free shaped tablet : heart, kidney, bean, water drop, half round, almond, flower shapes, etc.
- Hard gelatin capsule : capsule No. 00, 0, 1, 2, 3, 4, 5



\* Inspection applicability of respective products is subject to actual product testing and confirmation by Enclony.

## CAPACITY


 5mm Round : 400,000 tablets/hr

 10mm Oblong : 200,000 tablets/hr


 #5 Capsule : 150,000 capsules/hr

## DEFECT ANALYZATION

### Tablet



### Capsule



**ENCLONY**

**ENCLONY Co., Ltd.**

www.enclony.com

No.1901,243,Digital-ro, Guro-gu, Seoul, Korea 08382

E-mail : sales1@enclony.com Tel : +82 2 855 0095 Fax : +82 70 8650 4810

**ENCLONY GmbH**

Carl-Legien-Straße 15, 63073 Offenbach am Main, Germany  
(Sirius Business Park, Offenbach)

**DARRON**

 [info@darron.fr](mailto:info@darron.fr) 

 +33 (0)4 50 69 59 91